



Better Living Now, Inc. America's First Choice For Healthcare Products™
 185 Oser Ave., Hauppauge, NY 11788

Phone 631-348-0032 Toll-Free 1-800-854-5729
 Fax 631-348-7704 Toll-Free 1-800-654-7515

E-mail info@BetterLivingNow.com orders@BetterLivingNow.com
 Web site www.BetterLivingNow.com Education.BetterLivingNow.com

Infusion Set Comparison

Straight Infusion Set Comparison								
Company	Animas® (a Johnson & Johnson Company)	Unomedical	Roche® Insulin Delivery Systems, Inc.	Roche® Insulin Delivery Systems, Inc.	Medtronic MiniMed, Inc.	Asante	Medtronic Minimed, Inc.	Medtronic MiniMed, Inc.
Model	Inset™	Contact Detach	Accu-Chek® Ultraflex	Accu-Chek® Rapid-D	Quick-set®	Conset	Sure-T	mio™
Product Picture								
Connector Type	Luer-Lock	Luer-Lock	Luer-Lock	Luer-Lock	Luer-Lock, Paradigm	Luer lock, Snap	Luer-lock Paradigm	Luer-Lock, Paradigm
Cannula Length	6, 9 mm	6, 8, 10 mm	8, 10 mm	6, 8, 10 mm	6, 9 mm	6, 9	6, 8, 10 mm	6, 9 mm
Gauge	27 gauge intro needle, 0.68mm (width) cannula	29 gauge cannula	27 guage intro needle, 25 gauge cannula	28 guage cannula	27 gauge intro needle, 22 gauge cannula	27	29 gauge	27 gauge intro needle, 0.68mm (width) cannula
Tubing Length	23, 43 inches 60, 109 cm	24, 31, 43 inches 60, 80, 109 cm	24, 31, 43 inches 60, 80, 109 cm	24, 31, 43 inches 60, 80, 109 cm	23, 43 inches 60, 109 cm	23, 43	18, 23, 32 in.	18, 23, 32 inches 46, 60, 80 cm
Serter	all-in-one	None	Link Assist	None	Quick-Serter	included	None	all-in-one
Disconnect Type	at site	at site	at site	secondary location	at site	at site	at site	at site
Features	Small all-in-one inserter Various colors	Teflon cannula		Ultra thin, Stainless steel cannula	Transparent window to see point of insertion, audible click to confirm insertion		steel cannula	Small all-in-one inserter Various colors

The information provided here is based on product labeling at the time it was published obtained from current user guides, test strip package inserts and manufacturer web sites
 Trademarks are the property of their respective owners
 Marketing highlighted in blue are favorable

For More Information, Call Customer Service 1-800-854-5729 or The Insulin Pump Group 1-877-262-2182



Better Living Now, Inc. America's First Choice For Healthcare Products™
 185 Oser Ave., Hauppauge, NY 11788

Phone 631-348-0032 Toll-Free 1-800-854-5729
 Fax 631-348-7704 Toll-Free 1-800-654-7515

E-mail info@BetterLivingNow.com orders@BetterLivingNow.com
 Web site www.BetterLivingNow.com Education.BetterLivingNow.com

Infusion Set Comparison

Angled Infusion Set Comparison							
Company	Animas® (a Johnson & Johnson Company)	Animas® (a Johnson & Johnson Company)	Animas® (a Johnson & Johnson Company)	Asante	Roche® Insulin Delivery Systems, Inc.	Roche® Insulin Delivery Systems, Inc.	Medtronic MiniMed, Inc.
Model	Comfort™	Comfort™ Short	inset™ 30	Comfort	Tender	Tender "mini"	Silhouette®
Product Picture							
Connector Type	Luer-lock	Luer-lock	Luer-lock	Luer-lock, Snap	Luer-lock	Luer-lock	Luer-lock, Paradigm
Cannula Length	17mm	13mm	13mm	13, 17mm	17mm	13mm	13, 17mm
Gauge	27 gauge intro needle 25 gauge cannula	27 gauge intro needle 25 gauge cannula	27 gauge intro needle 0.68 (width) cannula		27 gauge intro needle 25 gauge cannula	27 gauge intro needle 25 gauge cannula	27 gauge intro needle 25 gauge cannula
Tubing Length	23, 31, 43 inches	23, 31, 43 inches	23, 43"	23, 43 inches	24, 31, 43 inches (60 cm, 80 cm, 110 cm)	24, 31, 43 inches (60 cm, 80 cm, 110 cm)	18, 23, 32, 43 inches
Serter	None	None	All-in-one	Manual	None	None	Sil-serter®
Disconnect Type	at site	at site	at site	at site	at site	at site	at site
Features			Small all-in-one inserter, various colors	Transparent window to see point of insertion, audible click to confirm insertion, angle insertion range is 20 to 45 degrees			

The information provided here is based on product labeling at the time it was published obtained from current user guides, test strip package inserts and manufacturer web sites
 Trademarks are the property of their respective owners
 Marketing highlighted in blue are favorable

For More Information, Call Customer Service 1-800-854-5729 or The Insulin Pump Group 1-877-262-2182