
# Prescribe With Confidence





## PARI LC SPRINT

The fastest addition to the PARI family of breath enhanced reusable nebulizers featuring consistently high respirable drug delivery PART #023F35



#### PARI LC PLUS

Breath enhanced reusable nebulizer considered the gold standard in nebulizer therapy for the past 10 years PART #022F81



### PARI LC STAR

Breath enhanced ultra fine particle nebulizer for effective penetration into smaller airways PART #022F51

PARI Reusable Nebulizers and compressors were used in pivotal drug trials by the following companies:



AstraZeneca DEY Genentech, Inc. **GILEAD INSPIRE Novartis SEPRACOR** 

PARI LC® reusable nebulizers feature breath enhanced technology resulting in increased aerosol delivery with decreased medication waste.







Performance Data of the PARI LC® Family of Reusable Nebulizers

	TOR	MMD	%RF	ATT
PARI LC® SPRINT	<b>500</b> mg/min	3.5 µm	68%	~ <b>5</b> min
PARI LC® PLUS	<b>440</b> mg/min	<b>3.8</b> µm	65%	<b>6-7</b> min
PARI LC® STAR	<b>400</b> mg/min	3.1 µm	<b>78</b> %	<b>6.5-8</b> min

TOR= TOTAL OUTPUT RATE [mg/min]

MMD= MASS MEDIAN DIAMETER [µm]

%RF= PERCENTAGE RESPIRABLE FRACTION [0.1~5 µm]

ATT= AVERAGE TREATMENT TIME [min]

\*Tested with the PARI PRONEB® Ultra II Compressor

<sup>1</sup>Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.

#### Breath-Enhanced Delivery

